NAHT DES MONATS: OKTOBER

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SUTURE OF THE MONTH: OCTOBER

OVARIOHYSTERECTOMY OF THE FEMALE DOG WITH Monosyn® & Novosyn®

The ovariohysterectomy (OHE) is a standard-surgery in the female dog in which ovaries and uterus are surgically removed. This reproductive surgery is often necessary to help control certain diseases influenced by reproductive hormones, e.g. mucometra, pyometra and cystic endometrial hyperplasia or to treat tumors, trauma or uterus torsion.

The requirements for a standard celiotomy are a dorsal recumbency of the anesthetized patient and a clipping and aseptically prepared ventral abdomen.

The incision should be performed directly caudal to the umbilicus. The length varies with the size of the dog. The incision goes through skin and subcutaneous tissue to expose the linea alba. After a careful stab incision through the linea alba the incision can be extended cranial and caudal to the stab with mayo scissors.

The elevation of the uterine horn can be performed with an ovariectomy hook or manually. Identify the suspensory ligament as a taut fibrous band at the cranial edge of the ovarian pedicle which fixates the ovaries caudal the kidneys. Stretch or break the suspensory ligament without tearing the ovarian vessels which are lying right underneath the band inside the fatty tissue. The ovarian artery has its origin in the aorta, so be sure to make a secure ligation. To perform the ligation use a clamp for fixation first.

The ideal suture materials for ligation are either absorbable, monofilament or polyfile (braided) sutures. With Monosyn® (monofilament) / or braided sutures you get a flexible suture material with easy handling and good knotting even in fatty tissues. The diameter depends on the size of the dog.

3: Settling of the ovary after ligation of the ovarian ligament

4: Check stump of short ovarian ligament for bleeding


The circumferential ligation with Novosyn® can be performed manually, using a Deschamps Needle. Alternatively a needle-thread combination (round needle) can be used.

We recommend a double circumferential ligation or a transfixing ligation (figure 8) or a combination of both.

After performing the ligatures cranial of the ovaries transect the ovarian pedicle distal to the clamp across the ovarian pedicle. Carefully remove the clamp and observe for hemorrhage. Separate the broad ligament from the uterine horns going caudally to the uterine body. Apply a Novosyn® ligation around the broad ligament if it is heavily infiltrated with vessels.


5: Placement of the ligature on the uterine body


6: Ligation uterine body completed, with clamp for protection

7: Overview cervical uterine ligation

1: Place the clamp in the cranial area of the uterine horn. A./V. ovarica must also be included

2: Placement of the ligature on the short ovarian ligament with A./V. ovarica

Ligate the uterine body cranial to the cervix. Depending on the diameter of the uterine body place a double circumferential ligation or transfixing ligation (figure 8) cranial the cervix. We recommend a separate ligation of the uterine vessels on both sides of the uterine body.

After placing a clamp across the uterine body cranial to the ligatures transect the uterine body cranial of the clamp. Remove the clamp from the uterine body and observe for hemorrhage like you have done at the ovary sides.

After replacing the uterine stump into the abdomen check for eventual hemorrhage. Close the abdominal wall in three layers (linea alba, subcutaneous tissue and skin) with a monofilament resorbable suture material. Braided suture material works like a wick and harbours the risk of bacterial invasion.

The recommendation for closing the abdominal wall is a monofilament long-term resorbable suture material which will be highlighted in Issue 3.

To close the subcutaneous tissue use Monosyn® with its atraumatic „tissue sliding“ in simple continuous suture patterns. These sutures are used to bring the tissue edges of the skin into anatomical apposition and to reduce tension.

You can also use Monosyn® for closing the skin. Even though this suture material is monofilament it has an excellent knotting quality so you can perform simple interrupted suture patterns or a Ford interlocking easily.

Monosyn® is also perfect for subcuticular closure. Good woundhealing also depends on choosing the right thread diameter and the applied tension to the tissue edges. So you can avoid ischemic islands of tissue and keep the inflammatory effects to a foreign body low.

(This text of performing a OHE is just a short essay of a complex surgery which can be performed in a variety of methods. Not all steps are named)

SOURCE:

Small Animal Surgery by Theresa Welch Fossum and own experience Dr. A. Wagner